


Gcomm MicroMotor



Product Categories: [ACCESSORIES](#)

Product Page: <http://flightdentalsystems.com/product/gcomm-micromotor/>

Product Summary

ST40 series made by G.Comm are brushed electric micro-motors, 10.000 - 40.000 rpm, with optional LED light and internal spray, detachable sleeve, compatible with Type E connection compliant with ISO 3964.

ST40 Micro-motors offer great performances with a smooth rotation up to 40.000 rpm and can be combined together with CMG01 controller unit, easily adaptable to all dental chairs equipped with a pneumatic or electric pedal.

Product Gallery



