

A6 Operatory Package



Product Categories: [OPERATORY SYSTEMS](#)

Product Tags: [a6 hybrid](#), [a6 operatory](#), [a6 operatory system](#), [a6 patient chair](#), [affordable](#), [brampton](#), [canada](#), [dental](#), [dental equipment](#), [electric dental chair](#), [flight dental](#), [flight dental systems](#), [hybrid dental chair](#), [hydraulic chair](#), [hydraulic dental chair](#)

Product Page: <https://flightdentalsystems.com/product/a6-operatory-package/>

Product Summary

The A6 Operatory System is a compact, ergonomic and modern dental operatory system. It features a hybrid lift system that combines a strong hydraulic base and a quiet electromechanical back motor that allows a smooth synchronized movement of the chair to deliver faster patient positioning.

In addition, the thin backrest and contoured base plate allows the user optimal maneuverability and access to the oral cavity. It comes standard with a one piece rotating cuspidor, Torch LED Light, Chair Controls for doctor and Assistant Side and a clean water bottle system.

Product Gallery







