

A6 Patient Chair



Product Categories: [PATIENT CHAIRS](#)

Product Page: <https://flightdentalsystems.com/product/a6-patient-chair/>

Product Summary

The A6 Patient Chair is a compact, ergonomic and modern dental chair. It features a hybrid lift system that combines a strong hydraulic base and a quiet electromechanical back motor that allows a smooth synchronized movement of the chair to deliver faster patient positioning. In addition, the thin backrest and contoured base plate allows the user optimal maneuverability and access to the oral cavity.

Product Gallery



