

A6 Patient Chair



Optimized by www.ImageOptimizer.net

Product Categories: [PATIENT CHAIRS](#)

Product Tags: [a6](#), [a6 operatory](#), [a6 patient chair](#), [adec](#), [affordable](#), [boyd](#), [brampton](#), [canada](#), [childrens chair](#), [dansereau](#), [dental](#), [dental chair](#), [dental equipment](#), [dental unit](#), [denturist chair](#), [electric chair](#), [electric dental unit](#), [electromechanical dental unit](#), [flight](#), [flight dental](#), [flight dental systems](#), [flight medical](#), [linak motors](#), [ortho chair](#), [pedo chair](#), [quiet dental chair](#), [quiet dental unit No](#)

Product Page: <https://flightdentalsystems.com/product/a6-patient-chair/>


Product Summary

The A6 Patient Chair is a compact, ergonomic and modern dental chair. It features a hybrid lift system that combines a strong hydraulic base and a quiet electromechanical back motor that allows a smooth synchronized movement of the chair to deliver faster patient positioning. In addition, the thin backrest and contoured base plate allows the user optimal maneuverability and access to the oral cavity.

Product Gallery



Optimized by [www.ImageOptimizer.net](#)



Optimized by [www.ImageOptimizer.net](#)