

FLIGHT CLAVE23



Product Categories: [AUTOCLAVES](#)

Product Tags: [9 Inch](#), [Autoclave](#), [canada](#), [Chamber](#), [Chamber Sterilizer](#), [chiropodist](#), [classN](#), [Clave16](#), [Clave23](#), [dental](#), [dental equipment](#), [Fast Sterilization](#), [flight](#), [flight dental systems](#), [Flight Sterilizer](#), [Health Canada Approved](#), [medical](#), [podiatry](#), [steam displacement](#), [Steam sterilizer](#), [sterilization](#), [tattoo](#), [veterinary](#)

Product Page: <https://flightdentalsystems.com/product/flight-clave23/>

Product Summary

The Clave23 Chamber Autoclave is a Class N Steam Sterilizer that uses state-of-the-art technology in steam sterilization using a pressure displacement technology that allows a more complete removal air inside the chamber before sterilization. In addition, this steam injection technology reduces maintenance, improves sterilization effectiveness and increases speed of the sterilization process.

The Clave23 uses the highest quality components and each unit is stress tested to ensure the quality and reliability of the machine. The Clave23 is compact and efficient sterilizer is suitable for any dental, veterinary, podiatry, chiropodist, tattoo, medical and beauty office.

As the name suggests this is the largest chamber autoclave offered by Flight Dental Systems at 23L in capacity.

Product Gallery



