

Mobile Cart - TRAD 2001



Product Categories: [DELIVERY SYSTEMS](#), [Uncategorized](#)

Product Tags: [dci](#), [delivery system](#), [dental](#), [dental chair](#), [dental equipment](#), [DENTAL FLIGHT](#), [dental handpieces](#), [dentist](#), [flight](#), [flight dental](#), [hygiene](#), [MC-1300F](#), [mobile cart](#), [mobile dental cart](#), [MOBILE DENTISTRY](#), [ORTHO CART](#), [pneumatic system](#), [SIDE CART](#), [TRAD 2001](#), [veterinary](#)

Product Page: <https://flightdentalsystems.com/product/mobile-cart-trad-2001/>

Product Summary

The Flight MC-1300F mobile cart is simple, economical and reliable. It has a U-Shaped Frame that provides sturdy and stable support for your delivery system, coupled with 4 dual casters that allows it to move smoothly on any surface. The MC-1300F Mobile cart has a large tray table for holding numerous instruments and cassettes adequately. The cart features an air driven pneumatic system that has individual air and water control for each handpiece as well as the syringe. There is also an extra slot for the optional Ultrasonic or Piezo Scaler to be added to this system. This cart is excellent for us in all Dental, Veterinary and Hygiene Operatories.

Product Gallery



