

NEW
COLLECTION

FLIGHT
DENTAL SYSTEMS

MK-dent[®]
GERMANY



ECO LINE

HIGH FLEXIBILITY, ECONOMIC PRICE!

- Smooth Grip
- Available with or without light
- Chrome coating
- Ceramic bearings
- Push Button system
- Glass rod light conductor

COUPLING-CONNECTION

HE21KL

High Speed Handpiece
with light

- 4-hole spray
- KaVo[®] Multiflex connection
- Standard Head
- 24W

HE22KL

High Speed Handpiece
with light

- 3-hole spray
- KaVo[®] Multiflex connection
- Small head
- 16 W

HE22K

High Speed Handpiece
without light

- 3-hole spray
- KaVo[®] Multiflex connection
- Small head
- 16 W

HB23KL

High Speed Handpiece
with light

- 1-hole spray
- KaVo[®] Multiflex connection
- Mini head
- 15 W

HB23K

High Speed Handpiece
without light

- 1-hole spray
- KaVo[®] Multiflex connection
- Mini head
- 15 W



Ceramic ball bearings
Longer durability compared to steel bearings



Glass rod light conductor
More than 40,000 LUX passing through



Orig. KaVo[®] compatible



Push button-System
Safe and easy to handle



NEW
Premium Spindle
Durability through perfect concentricity and carbide insert

**NEW
COLLECTION**

QUICK CONNECTORS



QC6016K

6-hole Quick Connector with LED

MK-dent & KaVo Multiflex®
connection



QC4014K

4-hole Quick Connector without light

MK-dent & KaVo Multiflex®
connection

AIR MOTORS



AM0004

4-hole Air motor without water supply

Max. 20,000 rpm

Standard Intramatic® connection,
suitable for all common instruments

ECO LINE LOW SPEED



Classic Grip

- Smooth Grip
- Available with or without light
- Chrome coating
- Intramatic® connection

LE11L | LE11

Contra Angle Handpiece
with light | without light

1:1, blue

With internal water

Max. 40,000 rpm

For CA burs(O 2.35 mm)



- KaVo® compatible
- Chrome coating
- Push button system

LCH11

Contra Angle Head

1:1, blue

Max. 40,000 rpm

For CA burs(O 2.35 mm)



Compatibility

MK-dent & KaVo Multiflex®
Low Speed Handpieces
and heads are compatible